



#### REFERENCES

- Drug master file No. 17416, FDA.
   Certification of Suitability of Monographs of the European Pharmacopoeia, EDQM.
   PRIME Journal Vol.9 Issue 2, 2019, P 74-76
   Aesthetic Plast Surg. 2017 Apr,41(2):422-428.

- A Aesthetic Plast Surg. 2017 Apr,41(2):422-428.
  Plast Reconstr Surg Glob Open 2020;8:e2648.
  Unpublished Clinical studies. Please request yvoire@lgchem.com
  Plast Reconstr Surg Glob Open. 2020 Jul 23;8(7):e2975.
  J Cosmet Dermatol. 2020 Dec;19(12):3229-3233.
  Data on file, LG Chem. Global shipment in 2020
  Data on file, LG Chem. Sales Data in 2022
  HA raw material of YVOIRE has equipped more strict quality requirement (Endotoxin, Protein level and sterility test results, etc) than European Pharmacopeia 9.0.
  US Food and Drug Administration. Executive Summary Dermal Filler Devices. November 18, 2008
  Carbohydr Polym. 2013;91(1):410-8.
  J Drugs Dermatol. 2011 Sep;10(9):974-80.
  Data on file, LG Chem.
  Galderma's patent, WO 2014/206701 A1
  ClinicalTrials.gov Identifier. NCT02334358
  ClinicalTrials.gov Identifier. NCT02661737
  ClinicalTrials.gov Identifier. NCT02106751

#### **(b)** LG Chem

E14, LG Science Park, 70, Magokjungang 10-ro, Gangseo-gu, Seoul 07795, Korea http://global.lgyvoire.com / yvoire@lgchem.com

These products are medical devices. Read the 'Precautions for Use' and 'How to Use' before use.

**LG** Chem

### LG aesthetics



## Basic Information of YVOIRE®

#### Product Summary

YVOIRE® provides diverse product options meeting the needs of doctors and patients.

Product		YVOIRE <sup>®</sup> classic	YVOIRE <sup>®</sup> volume	YV
		CLASSIC +	VOLUME +	CONTOUR
Practical uses		Smoothing and Filling Fine Wrinkles	Correcting deep wrinkles & Enhancing soft tissue volume	R &
Product image		VVOIRE CLASSIC	YVOIRE" VOLUME Wollwe Wollwe	VVOIRE CONTOUR CONTOUR CONTOUR CONTOUR CONTOUR
Composition	Relative particle size	Cross-linked Hyaluronic acid	Cross-linked Hyaluronic acid	Cros
	Lidocaine	0.3% lidocaine	0.3% lidocaine	-
Enclosed needles		30G x ½" 1ea 27G x ½" 1ea	27G x ½" 2ea	23G x 1" 1ea
Volume per syringe		1ml	1ml	1ml
Volumizing capacity				

### OIRE<sup>®</sup> CONTOUR



CONTOUR +

## **YVOIRE**<sup>®</sup> Product Information

YVOIRE<sup>®</sup> has natural effiacy due to delicate production.

HA raw material manufactured according to GMP-Rules and listed in US FDA DMF<sup>1</sup> with Certificate of suitability to the monographs. of the European Pharmacopeia by EDQM<sup>2</sup>

> FDA DMF COS

\* HA: Hyaluronic acid DMF: Drug Master File; COS: Certificates of Suitability

Improved efficacy and safety through HICE cross-linking technology<sup>3</sup>



\* HICE : High Concentration Equalized

Global sales in over 40 countries<sup>9</sup>



• UK, France, Germany, China, Thailand, Brazil, Russia and other countries • The countries shown on the map may not be accurate

One of few company that fully controls the quality and safety from HA raw material to end products



Confirmed effectiveness and safety through multiple clinical trials<sup>4,5,6,7,8</sup>



(more than 1,500 cumulated subjects)

More than 10 million YVOIRE injection globally.10



YVOIRE

### HA Features of **YVOIRE**<sup>®</sup>

History of LG Chem's HA research and development

HA raw material of YVOIRE<sup>®</sup> is directly manufactured by LG Chem and controlled with high standards as pharmaceutical products.<sup>11</sup> The HA raw material of YVOIRE<sup>®</sup> is quality assured, identical to the ophthalmic HA fine to put in eyes.





LG Chem's non-animal origin hyaluronic acid has been marketed in various medical fields since 1995 with confirmed safety and efficacy

LG Chem's Hyaluronic acid products

YVOIRE<sup>®</sup> is produced with HA raw material which has a high molecular weight of 3 million Daltons.<sup>3</sup>

<b>YVO</b> IRE <sup>°</sup>	3 million Daltons
Product A	1 million Daltons <sup>12</sup>
Product B	2.5 million Daltons <sup>12</sup>





## Cross-linking technology of YVOIRE<sup>®</sup>

#### HICE<sup>™</sup> Cross-linking technology of YVOIRE<sup>®</sup>

HICE<sup>™</sup> which stands for 'High Concentration Equalized' is a cross-linking technology developed by LG Chem. HICE<sup>™</sup> cross-linking technology facilitates a higher cross-linking ratio with few amount of cross-linkers by dispersing the agents into highly concentrated hyaluronic acid during crosslinking reaction.



#### The products chosen for the comparison had the highest storage modulus among its lidocaine options in the brand

a) Modification Degree: The MoD is calculated with the following formula, (Pendent cross-linkers + cross-linked cross-linkers) / (HA units) This represents the total cross-linkers(BDDE) in the end product only

b) Cross-linking Ratio: The number of cross-linked BDDE among the total BDDE

c) Storage modulus (G', Pa): Measured at 0.02Hz (Frequency sweep mode, 25°C, strain 4%), clinically indicates volume restoration effect<sup>14</sup>

d) Pa, Pascal, Storage modulus G'

#### Benefits of HICE<sup>™</sup> Cross-linking Technology of YVOIRE<sup>®</sup>



Cross-linking ratio of HA filler products<sup>16</sup>







### HA features of YVOIRE<sup>®</sup> CONTOUR

 $\mathsf{YVOIRE}^{\textcircled{B}}$  contour shows smooth surface outlines and supporting capacity owing to its ideal multi-particle size composition.  $^{3,5}$ 



YVOIRE<sup>®</sup> contour shows smooth gel surface by filling the void between main large particles with small particles, leads to high plasticity.



#### Microscopic view of YVOIRE<sup>®</sup>contour

Small particles fill the void between main large particles

Main large particles are the core of lifting capacity -





### Clinical Data of **YVOIRE**®

YVOIRE<sup>®</sup> confirmed its clinical efficacy on wrinkle correction and volume restoration, durability and safety through multiple clinical trial (more than 1,500 cumulative subjects).

#### 2008

A randomized, multi center, single-blind, activecontrolled, matched pairs design clinical study to evaluate the correction of facial wrinkles and folds and safety of Meso Glow versus IAL system in nasolabial fold intradermal injection

41 subjects, Korea

#### -----2009

A randomized, multi center, single-blind, activecontrolled, matched pairs design clinical study to evaluate the correction of facial wrinkles and folds up to 12 months after injection folds and safety of HA IDF versus Restylane<sup>®</sup> in nasolabial fold intradermal injection<sup>7</sup>

#### 58 subjects, Korea

#### 2010

A randomized, multi center, single-blind, activecontrolled, matched pairs design clinical study to evaluate the correction of facial wrinkles and folds and safety of HA IDF II versus Perlane<sup>®</sup> in nasolabial fold injection<sup>5</sup>

#### 57 subjects, Korea

#### 2012

A multi center, prospective and retrospective observation study to evaluate the safety and biodegradability of YVOIRE® classic s injected in the nasolabial folds up to 12 months after injection

.

68 subjects, Korea

#### 2013

Evaluation on safety and biodegradability of YVOIRE<sup>®</sup> volume injection in the nasolabial

80 subjects, Korea

#### 2013

A randomized, multi center, double-blind, active-controlled, matched pairs design clinical study to evaluate the efficacy and safety of HA IDF II plus versus HA IDF II in nasolabial fold injection<sup>8</sup>

62 subjects, Korea

Indicated year means the beginning year of clinical study, and the number of subjects is based on enrollment, Comparators are different in each clinical study, corresponds to each YVOIRE<sup>®</sup> products.

#### 2014

A randomized, multi center, double-blind, active- A multicenter, extension clinical study to evaluate the long-term efficacy and safety of controlled, matched pairs design clinical study repeat treatment of YVOIRE® contour injected to evaluatae the efficacy and safety of HA IDF plus versus HA IDF in nasolabial fold injection<sup>4</sup> into the anteromedial malar region in subjects who have completed the LG-HACL014 study

62 subjects, Korea

#### 2014

A multi-center, prospective and observational A multicenter, randomized, rater-blinded, activecontrolled, matched pairs design clinical study study to evaluate long-term safety and to evaluate the efficacy and safety of YVOIRE® biodegradability of modified sodium hyaluronate gel (YVOIRE<sup>®</sup> volume s) injected for correction contour compared with Restylane SubQ<sup>™</sup> injected into the anteromedial malar region<sup>5</sup> of nasolabial folds Medical Device Final Clinical Study Report<sup>18</sup>

83 subjects, Korea

#### 2014

A multi-center, prospective and observational study to evaluate long-term safety and biodegradability of YVOIRE® classic s injected into the nasolabial folds<sup>17</sup>

•

503 subjects, China

#### 2015

13 subjects, Korea \* Repeated Treatment

#### 2015

503 subjects, China

#### 2022

Additional global clinical trials are in progress

To Evaluate the Performance and Safety of YVOIRE<sup>®</sup> Volume Plus for Improvement of Midface Volume<sup>19</sup>

103 subjects, Germany

A Randomized, Multicenter, Evaluator-Blinded, Active-Controlled, Parallel-Group Investigation to Evaluate the Performance and Safety of YVOIRE<sup>®</sup> Classic Plus Versus Comparator for Temporary Correction of Nasolabial Folds<sup>20</sup>

104 subjects, Germany

## Clinical Data of **YVOIRE**<sup>®</sup>

#### Efficacy on wrinkle correction

YVOIRE<sup>®</sup> Classic Plus and Volume Plus received satisfying responses from overall subjects after 26 weeks from injection



\* Evaluation scale : Evaluate the degree of wrinkle correction by GAIS (global aesthetic improvement scale, -1~3 score) at 26 weeks after the injection \*\* above graph shows the number of subjects evaluated on each GAIS score (%)

#### Study design

#### Classic Plus

- A randomized, multicenter, double-blind, split-face study
- Total 62 subjects

• To evaluate the pain relief, efficacy of 26-weeks, and safety of 52-weeks after injecting YVOIRE<sup>®</sup> classic s and YVOIRE® classic plus for correction of nasolabial folds (50% linear threading technique, 50% serial puncture technique)<sup>4</sup>

#### Volume Plus

- A randomized, multicenter, double-blind, split-face study
- Total 62 subjects
- To evaluate the pain relief, efficacy of 26-weeks, and safety of 52-weeks after injecting YVOIRE® volume s and YVOIRE® volume plus for correction of nasolabial folds (100% linear threading technique)<sup>8</sup>

#### Pain relief

#### YVOIRE<sup>®</sup>Classic Plus and Volume Plus significantly reduced the patient's pain<sup>4,8</sup>





Subjects evaluated the pain using the VAS (Visual analogue scale) immediately after injection. The score ranges from 0-100 and a higher score indicates greater pain intensity.

### Clinical Data of **YVOIRE**<sup>®</sup>

#### Long-term Efficacy of YVOIRE<sup>®</sup>

YVOIRE<sup>®</sup> contour confirmed the clinical efficacy and safety of volume augmentation in mid-face for up to 1 year.<sup>5</sup>

#### Phase IV - YVOIRE<sup>®</sup> CONTOUR

Amulticenter, randomized, rater-blinded, active-controlled, Split-Face Design clinical study to evaluate the efficacy and safety of YVOIRE® contour (test device) compared with Product A" (Biphasic HA, control device) injected into the anteromedial malar region (n=83).



#### Long-term Safety of YVOIRE<sup>®</sup>

YVOIRE<sup>®</sup> classic s and YVOIRE<sup>®</sup> volume s confirmed the long term safety and biodegradability.<sup>17,18</sup>

#### Phase IV - YVOIRE<sup>®</sup> CLASSIC(S)

A multi-center, prospective and observational study to evaluate long-term safety and biodegradability of YVOIRE<sup>®</sup> classic s for 24 months after the initial injection into nasolabial folds.

- Target population : Chinese subjects whose WSRS grade is equal or greater than 2
- No. of subject : Total 503
- All the serious adverse events(1.99%) were unrelated to the investigational medical device.

The complete biodegradation was observed in all the subjects within 78 weeks after last treatment.

#### Phase IV - YVOIRE<sup>®</sup> VOLUME(S)

A multi-center, prospective and observational study to evaluate long-term safety and biodegradability of YVOIRE® volume s for 24 months after the initial injection into nasolabial fold.

- Target population : Chinese subjects whose WSRS grade is equal or greater than 2
- No. of subject : Total 503
- All of the serious adverse events (3.98%) are unrelated with the study device.

AE: adverse events

• The complete biodegradation was observed in all the subjects within 78 weeks after last treatment.

### Injection Area of **YVOIRE**<sup>®</sup>



These are for your information only, the injection areas and the appropriate products for it must be chosen according to the judgement of a trained physician. 'Indication' and 'Precautions for Use' may vary from countries.

\* Specific regions such as periobital and lips are not recommended for some countries(China,etc)







### Science for your Dream







# YVÖIRE®

### **Use Chem** LG aesthetics